


# Avalon Dual Lumen Cannula Dimensions



Description	Units	A	B			C		D		E		F		G		H		
		Connect or Size	Body Diameter			Insertable Length		Proximal Insertable Body to Infusion Port		SVC to Infusion Port Length		Infusion Port to Tip Length (Lead Length)		SVC to Tip Length		Lead Diameter		
		in	Fr.	mm	in	cm	in	cm	in	cm	in	cm	in	Fr.	in	Fr.	mm	in
Description Size	Order Code																	
13 Fr.	10013	1/4	13	4.3	0.17	11	4.3	8	3.2	1.7	0.7	2.8	1.1	4.5	1.8	11	3.7	0.15
16 Fr.	10016	1/4	16	5.3	0.21	14	5.5	10	3.9	3.1	1.2	4	1.6	7.1	2.8	14	4.7	0.19
19 Fr.	10019	1/4	19	6.4	0.25	21	8.3	15	6	4.3	1.7	5.7	2.3	10	4	16	5.3	0.21
20 Fr.	10020	3/8	20	6.6	0.26	31	12.2	22	8.5	5.3	2.1	9.4	3.7	14.7	5.8	17	5.7	0.22
23 Fr.	10023	3/8	23	7.7	0.3	31	12.2	22	8.5	5.3	2.1	9.4	3.7	14.7	5.8	20	6.7	0.26
27 Fr.	10027	3/8	27	9	0.35	31	12.2	22	8.5	5.3	2.1	9.4	3.7	14.7	5.8	24	8	0.31
31 Fr.	10031	3/8	31	10.3	0.41	31	12.2	22	8.5	5.3	2.1	9.4	3.7	14.7	5.8	27	9	0.35